
Next Generation

Pregnancy Testing IDEXX

dairyexpress.une.edu.au



Dairy Express offers a unique type of milk pregnancy testing through your herd testing samples. Original milk pregnancy testing was done through testing with milk progesterone.

The current generation of milk pregnancy testing is done by testing Pregnancy Associated Glycoprotein's (PAGS).

PAGS are a group of proteins produced by the maternal and embryonic regions of the placenta and therefore are specifically associated with pregnancy and are found in milk and serum.

This award winning laboratory test is manufactured by IDEXX and enables Dairy Express to detect PAGS in herd testing milk samples making it easy to detect pregnancy early and conduct routine pregnancy confirmations for your herd.



ECONOMIC BENEFITS & ADVANTAGES

- 35 days post breeding60 days post calving
- 98.8% sensitivity97.4% specificity
- Convenient & efficient
- Flexible
- ID integrity
- Reduces animal handling and stress
- Avoids lost production
- Helps reduces calving intervals

For additional information please contact Dairy Express

FREECALL 1800 633 165